

Product Datasheet

TNF ALPHA antibody (orb435375)

Description	Mouse monoclonal antibody to TNF ALPHA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, FA, IHC-Fr
Immunogen	Recombinant human tumour necrosis factor alpha produced in Escherichia coli.
Target	TNF ALPHA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.