

Product Datasheet

CD309 antibody (orb435376)

Description	Mouse monoclonal antibody to CD309
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant fragment corresponding to amino acids 1158 to 1345 within the intracellular region of mouse CD309.
Target	CD309
Preservatives	0.1% Sodium azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.