

Product Datasheet

ELMO1 antibody (orb435384)

Description	Goat polyclonal antibody to ELMO1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide with the sequence PKEPSNYDFVYDCN from the C terminus of the protein.Sequence according to the Entrez Gene accession number NP_055615; NP_56979
Target	ELMO1
Preservatives	0.02% Sodium Azide (NaN ₃)0.5% Bovine Serum Albumin. Tris buffered saline
Form/Appearance	Ig fraction - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.