

Product Datasheet

INFLUENZA HEMAGGLUTININ antibody (orb435413)

Description	Chicken polyclonal antibody to INFLUENZA HEMAGGLUTININ
Species/Host	Gallus
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	Synthetic peptide YPYDVPDYA from influenza hemagglutinin (HA).
Target	INFLUENZA HEMAGGLUTININ
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.