

Product Datasheet

KINESIN HEAVY CHAIN antibody (orb435417)

Description	Mouse monoclonal antibody to KINESIN HEAVY CHAIN
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, IHC-Fr, IHC-P, IP, WB
Immunogen	Bovine brain kinesin.
Target	KINESIN HEAVY CHAIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.