

## Product Datasheet

### CFTR antibody (orb435420)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CFTR                                                        |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IH-Fr, IP, WB                                                                            |
| <b>Immunogen</b>           | Fusion protein containing residues 1197-1480 of human CFTR ( <b>Chang et al. 1993</b> ). |
| <b>Target</b>              | CFTR                                                                                     |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> ). Phosphate Buffered Saline                         |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                    |
| <b>Concentration</b>       | 1.0 mg/ml                                                                                |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                                 |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | IgG1                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>Purity</b>              | Purified                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |