

Product Datasheet

HEPARIN/HEPARAN SULFATE antibody (orb435426)

Description	Mouse monoclonal antibody to HEPARIN/HEPARAN SULFATE
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA, RIA, WB
Target	HEPARIN/HEPARAN SULFATE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate Buffered Saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.