

Product Datasheet

INTEGRIN ALPHA 9 BETA 1 antibody (orb435429)

Description	Mouse monoclonal antibody to INTEGRIN ALPHA 9 BETA 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IP
Immunogen	Murine L cells transfected with human alpha-9 integrin.
Target	INTEGRIN ALPHA 9 BETA 1
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.