

Product Datasheet

DRP-3 antibody (orb435432)

Description	Rabbit polyclonal antibody to DRP-3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	Peptide with the sequence YDGPVFDLTTPK.
Target	DRP-3
Preservatives	0.1% Sodium Azide (NaN ₃)
Form/Appearance	Serum - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Serum
Expiration Date	12 months from date of receipt.