

Product Datasheet

GREEN FLUORESCENT PROTEIN antibody (orb435434)

Description	Mouse monoclonal antibody to GREEN FLUORESCENT PROTEIN
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Immunogen	GFP-BSA conjugate.
Target	GREEN FLUORESCENT PROTEIN
Preservatives	Bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.