

Product Datasheet

LIPID A LPS antibody (orb435440)

Description	Goat polyclonal antibody to LIPID A LPS
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Whole cell preparation of lipid A from E. coli O157.
Target	LIPID A LPS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	4.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.