

Product Datasheet

POLIOVIRUS antibody (Biotin) (orb435442)

Description	Goat polyclonal antibody against POLIOVIRUS conjugated to Biotin
Species/Host	Goat
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF
Immunogen	Pool of purified poliovirus (types 1-3)
Target	POLIOVIRUS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	IgG fraction conjugated to Biotin - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.