

Product Datasheet

TRYPANOSOMA CRUZI LPG antibody (orb435453)

Description	Mouse monoclonal antibody to TRYPANOSOMA CRUZI LPG
Species/Host	Mouse
Reactivity	Protozoa
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, ICC, IF
Immunogen	T. cruzi LPG purified from epimastigotes.
Target	TRYPANOSOMA CRUZI LPG
Preservatives	None present
Form/Appearance	Ascitic fluid - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Ascites
Expiration Date	12 months from date of receipt.