

Product Datasheet

INDOLEAMINE 2,3-DIOXYGENASE antibody (orb435458)

Description	Mouse monoclonal antibody to INDOLEAMINE 2,3-DIOXYGENASE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	Peptide corresponding to amino acids 78-184 of human IDO fused to GST.
Target	INDOLEAMINE 2,3-DIOXYGENASE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.