

Product Datasheet

GLUCOSE TRANSPORTER 1 antibody (orb435463)

Description	Rabbit polyclonal antibody to GLUCOSE TRANSPORTER 1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	12 amino acid peptide of the C-terminus of hGLUT-1 (LFHPLGADSQV) conjugated to KLH.
Target	GLUCOSE TRANSPORTER 1
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.