

Product Datasheet

Seroblock Fcr antibody (orb435575)

Description	Mouse monoclonal antibody to Seroblock Fcr
Species/Host	Mouse
Reactivity	Mouse
Tested Applications	Blocking
Immunogen	PU5 1.8 IOE7 Balb/c mouse cell line.
Target	MOUSE SEROBLOCK FCR
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Expiration Date	12 months from date of receipt.