

Product Datasheet

Bacterial LIPOTEICHOIC ACID protein (orb435584)

Description	Native bacterial LIPOTEICHOIC ACID protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	LIPOTEICHOIC ACID
Preservatives	None present.
Form/Appearance	Purified lipoteichoic acid - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.