

Product Datasheet

Viral HEPATITIS B 320 CORE ANTIGEN protein (orb435587)

Description	Recombinant of viral HEPATITIS B 320 CORE ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	HEPATITIS B 320 CORE ANTIGEN
Preservatives	50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Application notes	18kDa AYW
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.