

Product Datasheet

Human ALPHA 1 ANTITRYPSIN protein (orb435591)

Description	Native human ALPHA 1 ANTITRYPSIN protein
Reactivity	Human
Tested Applications	ELISA
Target	ALPHA 1 ANTITRYPSIN
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified alpha 1 antitrypsin from human serum - liquid
Concentration	2 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.