

Product Datasheet

Human COLLAGEN IV protein (orb435600)

Description	Native human COLLAGEN IV protein
Reactivity	Human
Tested Applications	ELISA
Target	COLLAGEN IV
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified protein - liquid
Concentration	1.0mg/ml.
Storage	+4°C, -20°C if preferred
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.