

Product Datasheet

Human FIBRIN DEGRADATION PRODUCT E protein (orb435618)

Description	Native human FIBRIN DEGRADATION PRODUCT E protein
Reactivity	Human
Tested Applications	ELISA
Target	FIBRIN DEGRADATION PRODUCT E
Preservatives	None present. Glycine buffered saline
Form/Appearance	Purified protein from human serum - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.