

Product Datasheet

Bacterial BETA GALACTOSIDASE protein (orb435621)

Description	Native bacterial BETA GALACTOSIDASE protein
Reactivity	Bacteria
Tested Applications	FA
Target	BETA GALACTOSIDASE
Preservatives	Contains Sucrose. Phosphate buffered saline
Form/Appearance	Purified protein from Escherichia coli - lyophilized
Concentration	Protein content is approximately 20-25% of total weight
Storage	-20°C after reconstitution
Note	For research use only
Purity	Lyophilized
Expiration Date	6 months from date of receipt.