

Product Datasheet

Bacterial MYCOPLASMA PNEUMONIAE protein (orb435639)

Description	Native bacterial MYCOPLASMA PNEUMONIAE protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	MYCOPLASMA PNEUMONIAE
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Inactivated Mycoplasma pneumoniae - liquid
Concentration	Current, batch-specific concentration 2.6 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.