

Product Datasheet

Bacterial PROTEIN G protein (orb435646)

Description	Recombinant of bacterial PROTEIN G protein
Reactivity	Bacteria
Tested Applications	ELISA, IP
Target	PROTEIN G
Preservatives	Recombinant Protein - lyophilized
Form/Appearance	Recombinant Protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.