

Product Datasheet

Human RETINOL BINDING PROTEIN protein (orb435656)

Description	Native human RETINOL BINDING PROTEIN protein
Reactivity	Human
Tested Applications	ELISA
Target	RETINOL BINDING PROTEIN
Preservatives	None present
Form/Appearance	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.