

Product Datasheet

Bacterial STREPTAVIDIN protein (orb435658)

Description	Bacterial STREPTAVIDIN protein
Species/Host	Bacteria
Tested Applications	ELISA
Target	STREPTAVIDIN
Preservatives	1% Mannitol1% Dextran. 0.01M Sodium Phosphate0.01M Sodium Borate0.15M Sodium Chloride
Form/Appearance	Purified streptavidin from Streptomyces avidinii - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Expiration Date	6 months from date of receipt.