

Product Datasheet

Human TROPONIN I protein (orb435665)

Description	Native human TROPONIN I protein
Reactivity	Human
Tested Applications	ELISA
Target	TROPONIN I
Preservatives	15mM Mercaptoethanol. 7M Urea5mM EDTA20mM Tris
Form/Appearance	Purified protein from human heart - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Application notes	(CARDIAC)
Purity	Purified
Expiration Date	6 months from date of receipt.