

Product Datasheet

Human TROPONIN T protein (orb435666)

Description	Native human TROPONIN T protein
Reactivity	Human
Tested Applications	ELISA
Target	TROPONIN T
Preservatives	None present. Urea/Tris buffer, pH 7.5 (7 M urea, 5 mM EDTA, 15 mM mercaptoethanol, 20 mM Tris)
Form/Appearance	Purified protein from human heart - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.