

Product Datasheet

Bacterial NEISSERIA GONORRHOEAE protein (orb435674)

Description	Bacterial NEISSERIA GONORRHOEAE protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	NEISSERIA GONORRHOEAE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Inactivated, homogenised Neisseria gonorrhoeae - frozen
Concentration	0.41 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.