

Product Datasheet

Bacterial CLOSTRIDIUM DIFFICILE TOXOID A protein (orb435678)

Description	Bacterial CLOSTRIDIUM DIFFICILE TOXOID A protein
Reactivity	Bacteria
Tested Applications	ELISA
Target	CLOSTRIDIUM DIFFICILE TOXOID A
Preservatives	None present. 0.05M HEPES, 0.15M NaCl, 5% Sucrose
Form/Appearance	Inactivated toxoid derived from purified Clostridium difficile Toxin A - lyophilized
Concentration	0.2 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.