

Product Datasheet

Viral EPSTEIN-BARR VIRUS NUCLEAR ANTIGEN protein (orb435680)

Description	Recombinant of viral EPSTEIN-BARR VIRUS NUCLEAR ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	EPSTEIN-BARR VIRUS NUCLEAR ANTIGEN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Recombinant Protein - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.