

Product Datasheet

Viral HEPATITIS B CORE ANTIGEN protein (orb435681)

Description	Recombinant of viral HEPATITIS B CORE ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA, WB
Target	HEPATITIS B CORE ANTIGEN
Preservatives	1mM EDTA. TRIS buffered saline, EDTA
Form/Appearance	Recombinant Protein - liquid
Concentration	0.1 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.