

Product Datasheet

Viral HEPATITIS B CORE ANTIGEN protein (orb435682)

Description	Recombinant of viral HEPATITIS B CORE ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	HEPATITIS B CORE ANTIGEN
Preservatives	50% Glycerol. Phosphate buffered saline
Form/Appearance	Recombinant protein produced in Escherichia coli - liquid
Concentration	1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.