

Product Datasheet

Viral HEPATITIS B SURFACE ANTIGEN AY protein (orb435683)

Description	Recombinant of viral HEPATITIS B SURFACE ANTIGEN AY protein
Reactivity	Virus
Tested Applications	ELISA, WB
Target	HEPATITIS B SURFACE ANTIGEN AY
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Recombinant Protein - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.