

Product Datasheet

Viral HEPATITIS B CORE ANTIGEN protein (orb435684)

Description	Recombinant of viral HEPATITIS B CORE ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA, WB
Target	HEPATITIS B CORE ANTIGEN
Preservatives	None present. Acetate buffer, pH 4.5
Form/Appearance	Recombinant protein produced in Pichia pastoris - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.