

Product Datasheet

Viral HEPATITIS C CORE ANTIGEN protein (orb435685)

Description	Recombinant of viral HEPATITIS C CORE ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	HEPATITIS C CORE ANTIGEN
Preservatives	8M Urea; 10mM betamercaptoethanol. 20mM Tris
Form/Appearance	Purified recombinant protein - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.