

Product Datasheet

Human ERYTHROPOIETIN ALPHA protein (orb435689)

Description	Recombinant of human ERYTHROPOIETIN ALPHA protein
Reactivity	Human
Tested Applications	FA
Target	ERYTHROPOIETIN ALPHA
Preservatives	None present.. Citrate buffered saline
Form/Appearance	Recombinant protein - lyophilized
Concentration	0.2mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.