

Product Datasheet

Human PROLACTIN protein (orb435692)

Description	Recombinant of human PROLACTIN protein
Reactivity	Human
Tested Applications	FA
Target	PROLACTIN
Preservatives	See Buffer solution.. Phosphate buffered saline
Form/Appearance	Recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.