

Product Datasheet

Protozoan PLASMODIUM VIVAX CSP PROTEIN protein (orb435699)

Description	Recombinant of protozoan PLASMODIUM VIVAX CSP PROTEIN protein
Reactivity	Protozoa
Tested Applications	ELISA
Target	PLASMODIUM VIVAX CSP PROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Recombinant Protein - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.