

Product Datasheet

Mouse ST6GALNAC2 protein (orb4357)

Description	Recombinant of mouse ST6GALNAC2 protein (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mMTris,150mMNaCl,pH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mMTris-HCl, 150mM NaCl, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Ser28-Arg373
Source	Human Cells
Uniprot ID	P70277
Expiration Date	6 months from date of receipt.