

Product Datasheet

Protozoan TOXOPLASMA GONDII P29 GRA7 protein (orb435701)

Description	Recombinant of protozoan TOXOPLASMA GONDII P29 GRA7 protein
Reactivity	Protozoa
Tested Applications	ELISA, WB
Target	TOXOPLASMA GONDII P29 GRA7
Preservatives	50% Glycerol. 50mM Tris-HCl1.5M Urea
Form/Appearance	Purified Recombinant Protein - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	(aa24-100)
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.