

Product Datasheet

Protozoan TOXOPLASMA GONDII P30 SAG1 protein (orb435702)

Description	Recombinant of protozoan TOXOPLASMA GONDII P30 SAG1 protein
Reactivity	Protozoa
Tested Applications	ELISA, WB
Target	TOXOPLASMA GONDII P30 SAG1
Preservatives	1.5 M Urea50% glycerol. 50 mM TRIS-HCl
Form/Appearance	Recombinant Protein - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	(aa45-198)
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.