

Product Datasheet

Chicken IL1 BETA protein (orb435708)

Description	Recombinant of chicken IL1 BETA protein
Reactivity	Gallus
Tested Applications	FA
Target	IL-1 BETA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified His-tagged recombinant protein expressed in E.coli - liquid
Concentration	250 µg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.