

Product Datasheet

Chicken IFN ALPHA protein (orb435710)

Description	Recombinant of chicken IFN ALPHA protein
Reactivity	Gallus
Tested Applications	FA
Target	IFN ALPHA
Preservatives	None present. 0.05M Sodium Acetate1M Sodium Chloride0.01M Tris-HCl2M Urea10mM Beta-Mercaptoethanol
Form/Appearance	Purified recombinant protein expressed in E.coli - liquid.
Concentration	1.0mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.