

Product Datasheet

Bovine IL4 protein (orb435716)

Description	Recombinant of bovine IL4 protein
Reactivity	Bovine
Tested Applications	ELISA, FA
Target	IL-4
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.04mg/ml
Storage	+4°C, -70°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.