

Product Datasheet

Bovine IL1 BETA protein (orb435718)

Description	Recombinant of bovine IL1 BETA protein
Reactivity	Bovine
Tested Applications	ELISA, FA
Target	IL-1 BETA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified Recombinant Protein - lyophilized
Concentration	40µg/ml
Storage	+4°C, -70°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Purified
Expiration Date	6 months from date of receipt.