

## Product Datasheet

### Bovine IL1 ALPHA protein (orb435736)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Recombinant of bovine IL1 ALPHA protein    |
| <b>Reactivity</b>          | Bovine                                     |
| <b>Tested Applications</b> | ELISA, FA                                  |
| <b>Target</b>              | IL-1 ALPHA                                 |
| <b>Preservatives</b>       | 10% Trehalose. Phosphate buffered saline   |
| <b>Form/Appearance</b>     | Purified recombinant protein - lyophilized |
| <b>Storage</b>             | -20°C only (ship +4°C)                     |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Recombinant                                |
| <b>Purity</b>              | Rec. Protein                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.             |