

Product Datasheet

Bovine LEUKEMIA INHIBITORY FACTOR protein (orb435741)

Description	Recombinant of bovine LEUKEMIA INHIBITORY FACTOR protein
Reactivity	Bovine
Tested Applications	ELISA, WB
Target	LEUKEMIA INHIBITORY FACTOR
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.