

## Product Datasheet

### Dog IgG protein (orb435745)

<b>Description</b>	Dog IgG protein
<b>Reactivity</b>	Canine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	IgG
<b>Preservatives</b>	Purified IgG - lyophilized
<b>Form/Appearance</b>	Purified IgG - lyophilized
<b>Concentration</b>	10.0 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	Purified
<b>Expiration Date</b>	6 months from date of receipt.