

Product Datasheet

Dog C-REACTIVE PROTEIN protein (orb435746)

Description	Recombinant of dog C-REACTIVE PROTEIN protein
Reactivity	Canine
Tested Applications	ELISA
Target	C-REACTIVE PROTEIN
Preservatives	None present. TRIS buffered saline, 2 mM CaCl ₂
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.